

EXOTICS MOSAIC

A *Saccharomyces cerevisiae* x *Saccharomyces paradoxus* hybrid.

A novel yeast for EXOTIC, COMPLEX, barrel fermented ICONIC WHITE WINES.

EXOTICS
MOSAIC

S. cerevisiae x *S. paradoxus* hybrid

ICONIC, COMPLEX AND
INTENSE WINES

ORIGIN

Exotics Mosaic is a product of the yeast hybridisation program of the Institute for Wine Biotechnology, Stellenbosch University, South Africa.

APPLICATION

Exotics Mosaic is a hybrid between *S. cerevisiae* and *S. paradoxus*. *S. paradoxus* is the closest relative to *S. cerevisiae* and is naturally found on grapes. This hybrid inherited the aromatic capabilities of both its parents, thereby expanding the aromatic potential and complexity of what pure *S. cerevisiae* have to offer. White wines produced by this yeast are described as having exotic aromas and flavours. Exotics Mosaic is recommended for the production of barrel fermented Chardonnay, Chenin blanc and Viognier. It enhances guava, passion fruit, tropical fruit salad and stone fruit aromas and flavours. It is cold sensitive and ferments at a steady rate in barrels, producing good mouth feel.

FERMENTATION KINETICS

Conversion factor: 0.56 - 0.61.

TECHNICAL PARAMETERS

Cold tolerance	18°C (64°F)
Optimum temperature range	18 - 20°C (64 - 68°F)
Osmotolerance	25 °Brix, 13.9 Baumè
Alcohol tolerance at 20°C (86°F)	15.5%

METABOLIC CHARACTERISTICS

Glycerol production	generally 9 -10 g/L
Volatile acidity production	generally lower than 0.4 g/L
SO ₂ production	none to very low
Nitrogen requirement	average
Malic acid degradation	none observed in white wines
Pectinolytic activity	yes

PHENOTYPE

Killer positive.

DOSAGE

30 g/hL (2.5 lb/1000 gal)

PACKAGING

Exotics Mosaic is vacuum-packed in 250 g packets. It must be stored in a cool (5 - 15°C, 41 - 59°F), dry place, sealed in its original packaging.



Anchor
OENOLOGY

Dedicated to fermentation excellence

DISTRIBUTOR:

OENOBRANDS SAS

Parc Agropolis II - Bât. 5
2196 boulevard de la Lironde
CS 34603 F-34397 Montpellier Cedex 5

info@oenobrands.com
www.oenobrands.com

RCS Montpellier - SIREN 521 285 304